

# VXS-115 VETERINARY X-RAY


**sri**  
source-ray, inc.  
REIMAGINING X-RAY




The VXS-115 Veterinary Digital X-Ray system is a high quality solution for reliable and convenient veterinary imaging. With flexible mobility, easy to use digital X-ray software, and compatibility for equine and bovine applications, this system is designed to meet all your veterinary needs.

## KEY ADVANTAGES

- ✓ Plugs into a standard wall outlet
- ✓ Easy to move around to clean or relocate
- ✓ Lubrication-free sliding system
- ✓ Next-Gen wired and wireless DR detectors
- ✓ Simple, easy to use digital X-ray software
- ✓ Limited 03 Year Warranty



**Compatible with Wired /  
Wireless X-ray Detector**



**User Friendly  
X-ray Acquisition Software**





## Model VXS-115 Specifications

### Generator

High frequency resonant inverter 1.5 kW

<b>Input Voltage:</b>	Input Voltage: 115/230 VAC, 50/60Hz
<b>kVp Range:</b>	40 - 100 kVp, adjustable in 1 KV increments closed loop regulation
<b>mA Range:</b>	15 Fixed. Constant mA, Independent of kV or time
<b>mAs Range:</b>	0.15 - 60
<b>Indicators:</b>	Digital display of kVp and mAs. X-Ray, Ready & Fault indicators.

### Exposure Switch

Detachable two-position

### X-Ray Tube

Stationary Anode 100 kVp

<b>Filtration:</b>	0.8 mm Al equivalent at 80 kV
<b>Target Material:</b>	Tungsten
<b>Anode Capacity:</b>	20,000 HU Storage
<b>Focal Spot:</b>	1.0 mm
<b>Target Angle:</b>	15°

### Collimator

Ralco R 72 S DHHS

<b>Laser Crosshairs:</b>	1.0 mm
<b>Tape Measure:</b>	SID Measurement
<b>Inclinometer:</b>	Angle Measurement

### Footswitch

Optional footswitch assembly for taking hands free exposures



Can be moved up, down or sideways along  
with the panel



MDSAP / ISO 13485 Certified

Specifications are subject to change without notice.

Names of companies and/or products appearing in this document are trademarks and/or registered trademarks of the respective companies.

Source-Ray, Inc.  
50 Fleetwood Court, Ronkonkoma, NY 11779  
(631) 244-8200  
x\_raysales@sourceray.com

Brochure, VXS-115, Rev.0